

Meerkat

User information appropriate for use by people with visual, reading or cognitive disability	The latest version of all instructions are always available on the R82 website support.R82.org and can be printed in larger sizes
Information available in pictogram	The instructions are made as pictograms
Intended use and intended environment including standard options	<p>The product is a standing frame for both anterior and posterior up right standing. The support given by using different accessories will compensate for the lack of function the child might have. As an accessory the Rocker can be mounted under the frame to make the product dynamic.</p> <p>The target group for this device is disabled children with cerebral palsy (CP). The consequence of CP is either spasticity, reduced walking or standing skills and/or impairments in addition to their motor disability. Other effects with no relation to this device is intellectual, vision and hearing impairment.</p> <p>Also children with muscular dystrophy, spina bifida, spinal cord injuries and neurological impairments can be subject to the device. An acknowledged classification of disabled persons, is the Gross Motor Function Classification System. The target group for this device is level II - IV.</p> <p>The product comes in 3 sizes. The product is an indoor product and can in the case of an emergency be taken outside.</p>

Maintenance instructions	<p>Maintenance of a medical device is wholly the responsibility of the owner of that device. Failure to maintain a device in accordance with the instructions may invalidate the device's warranty. Furthermore, failure to maintain a device may compromise the clinical condition or safety of users and/or their caregivers.</p>											
<table border="0"> <tr> <td>Frequency</td> <td>Description</td> </tr> <tr> <td>Before every use</td> <td>- All nuts and bolt on this product must be checked and tightened to avoid failures</td> </tr> <tr> <td>Daily</td> <td>- Use a dry cloth to clean the product</td> </tr> </table>	Frequency	Description	Before every use	- All nuts and bolt on this product must be checked and tightened to avoid failures	Daily	- Use a dry cloth to clean the product	<table border="0"> <tr> <td>- Visual inspection. Repair or replace any damaged or worn-out parts</td> </tr> <tr> <td>- Ensure that all fasteners are in place and fitted correctly</td> </tr> <tr> <td>- Check all fittings and buckles for signs of worn out parts</td> </tr> <tr> <td>- Be aware of any MAX mark indications before adjustment</td> </tr> <tr> <td>- Wipe off any residue and dirt from the product using a cloth with warm water and a mild detergent/soap and dry before use</td> </tr> </table>	- Visual inspection. Repair or replace any damaged or worn-out parts	- Ensure that all fasteners are in place and fitted correctly	- Check all fittings and buckles for signs of worn out parts	- Be aware of any MAX mark indications before adjustment	- Wipe off any residue and dirt from the product using a cloth with warm water and a mild detergent/soap and dry before use
Frequency	Description											
Before every use	- All nuts and bolt on this product must be checked and tightened to avoid failures											
Daily	- Use a dry cloth to clean the product											
- Visual inspection. Repair or replace any damaged or worn-out parts												
- Ensure that all fasteners are in place and fitted correctly												
- Check all fittings and buckles for signs of worn out parts												
- Be aware of any MAX mark indications before adjustment												
- Wipe off any residue and dirt from the product using a cloth with warm water and a mild detergent/soap and dry before use												
Weekly	<table border="0"> <tr> <td>- Check that all castors are moving freely and that all wheel locks/brakes are easily activated</td> </tr> <tr> <td>- Disinfect the product. please check the disinfection section in this User Guide for more information</td> </tr> </table>	- Check that all castors are moving freely and that all wheel locks/brakes are easily activated	- Disinfect the product. please check the disinfection section in this User Guide for more information									
- Check that all castors are moving freely and that all wheel locks/brakes are easily activated												
- Disinfect the product. please check the disinfection section in this User Guide for more information												
Every 6 months	<table border="0"> <tr> <td>- Check the condition of the upholstery and replace it if necessary</td> </tr> </table>	- Check the condition of the upholstery and replace it if necessary										
- Check the condition of the upholstery and replace it if necessary												
Annually	<table border="0"> <tr> <td>- Inspect the frame for cracks or signs for worn-out parts and carry out annual service. Never use a product which appears to be substandard or faulty</td> </tr> <tr> <td>- The product has to be checked and tightened by a qualified technician at least every 12 months to avoid failures.</td> </tr> </table>	- Inspect the frame for cracks or signs for worn-out parts and carry out annual service. Never use a product which appears to be substandard or faulty	- The product has to be checked and tightened by a qualified technician at least every 12 months to avoid failures.									
- Inspect the frame for cracks or signs for worn-out parts and carry out annual service. Never use a product which appears to be substandard or faulty												
- The product has to be checked and tightened by a qualified technician at least every 12 months to avoid failures.												

Cleaning description

Washing:

Meerkat

The cushioning is removable and the cover can be washed at max. 60° using a mild detergent in a washing machine.

Warewashing machine

This product can be washed at 60° using a mild detergent for 15-20 minutes in a warewashing machine designed for medical devices. Use the machine's drying function to dry the product. Remove all accessories before washing the main product. Clean the accessories separately. Please allow the product to dry before reuse.

Disinfection description

The product can be disinfected with a 70% IPA solution. It is recommended to wipe off any residue and dirt from the product, using a cloth wrung in warm water and a mild detergent/soap without chlorine and let it dry before disinfection.

Dimensions

<http://www.R82.org> - Select product

Mass

<http://www.R82.org> - Select product

Maximum slope

<http://www.R82.org> - Select product

Type(s) of tyres

Castors

Used in combination with other products

N/A

Warning about dangerous combinations

The CE mark must be removed when rebuilding the product, when used in combination with a product from another manufacturer or when using other than original R82 spare parts and fittings.

The Meerkat might lose its flame resistant characteristics if using other than original R82 cushioning.

List of accessories, detachable parts and materials

<http://www.R82.org> - Select product

Programmable controller information

N/A

Operator control adjustments

Manual

Folding and dismantling

See pictograms in the User Guide

Used in motor vehicles

N/A

Transport instructions (unoccupied)

Land and/or air

Sound power level

N/A