

## HEALTH CORNER

VOLUME 3, ISSUE 1 - December 2019

## Do you have a client who would benefit from the Salute RDX?

Dynamic Alternating Pressure
Therapy creates an area of zero
pressure leading to
"True Tissue Off Load",
inducing hypermeia- a critical
component to increasing
healing time

No "bottom out" meaning it will remain inflated even in a power outage

> The "Happy Heel Protection" system bleeds out 25% of air from the heel section in order to soften up the heel area for fragile heels



Can be a cost effective solution over several other mattress types

No periodic checking and/or adjustments for inflation

## **Salute RDX**

HME's team is a class leader in surface and low air loss mattress Provides low air loss,
which assists microclimate
managament
(temperature and moisture
control) to keep the wound
dry, cool, and promote healing

Victoria Location 2521 Government Street Victoria, BC, V8T 4P6 Office Phone: 250-386-0075 Richmond Location 130-4011 Viking Way Richmond, BC, V6V 2K9 Office Phone: 604-821-0075 Surrey Location 140-19288 22 Ave Surrey, BC, V3Z 3S6 Office Phone: 604-535-5768

Education, Leadership, Health Promotion & Innovation

www.hmebc.com